


ClickTech Oligo Link Kit - Dye Labeling Reactions of Oligos

This protocol was especially tested for performing click chemistry with 70 pmol to 90 nmol of modified oligonucleotide using our proprietary solid catalyst (the “reactor”) (available as single product Solid Reactor BCMI-008). As the reactor won’t be dissolved during the reaction, the reaction set up is extremely easy, the reaction faster and purification simplified. Please be aware, that this protocol is only meant as a starting point. For other amounts and reaction partners used, please consider the user manual of the ClickTech Oligo Link Kits (BCK-OL) to obtain the optimal reaction outcome.



You will need following reagents and equipment:

- Alkyne- or azide-modified oligonucleotide
- 10 mM solution of your preferred label-azide/ alkyne (dye, biotin, targeting moieties and many more)
- Eventually DMSO/H₂O to dissolve your label-azide in case you bought is in solid form
- ClickTech Oligo Link kit (BCK-OL)
- Table centrifuge
- Thermomixer
- Purification (e.g. ethanol precipitation, HPLC, purification kits...)
- Analytical HPLC system

Considerations:

- The “Reactor” contains a stable **heterogeneous catalyst**, which won’t be dissolved during the reaction.
- The click reaction can be performed with 10-100 μ M DNA oligonucleotide solutions using this basic click protocol. For more concentrated samples a “preparative click” protocol might be needed. For RNA oligonucleotides check extra section in the user manual.
- Only terminal alkynes can react with azides using the kit reaction conditions.

Click reaction procedure:

1. Dissolve/Dilute your label azide/label alkyne in DMSO or water to 0.2-10 mM to use it as a stock solution.
2. Add the appropriate amount of 10x Activator² to the Reactor, e.g. 2.5 μL 10x Activator² are added to Reactor 25 to be used at a total reaction volume of 25 μL . Depending on Reactor amount and final volume, this needs to be adjusted (see Table 1).
3. Add the alkyne/azide modified DNA oligonucleotide to the vial to a final concentration of 10-100 μM .
4. Add 2 equivalents of label azide/alkyne per equivalent of alkyne/azide in the oligonucleotide. For example, a 10 μM solution of a singly alkyne-modified oligonucleotide is mixed with 20 μM of a label azide for the click reaction.
5. Close the vial and incubate the mixture at 45 $^{\circ}\text{C}$, 600 rpm for 1 h in a thermomixer. Alternatively, a water bath can be used. When using fluorophores, protect the vial from light. Make sure that the Reactor is within the reaction solution during the reaction. Spin down the solution if needed.
6. After the reaction: Spin down the Reactor. Transfer the supernatant with the clicked oligonucleotide to a new vial. **Note:** For long-term storage, reacted samples (without Reactor) should be kept at -20°C .
7. Analyze the reaction mixture by gel electrophoresis, HPLC or ion exchange chromatography (IEC). Purifications using column-based kits for oligonucleotide purification (e.g. PCR purification kit from Qiagen) give good results. Make sure the length of the oligonucleotide is compatible with the purification kit. Alternatively, purification can be done by HPLC, IEC or ethanol precipitation (product loss likely).

Exemplary Label-Oligo Click

This guide will help you decide which stock solution concentration of the label azide/alkyne should be prepared. All concentrations within the table refer to the stock solution concentrations for exemplary setups using two equivalents of label azide/alkyne for a singly modified oligonucleotide alkyne/azide.

Table 1: Exemplary volumes needed for reaction setups of “basic” label-oligonucleotide click reactions.

Reactor	V (Activator ²)	c (Oligo)	V (Oligo)	c (Label)	V (Label)	V (H ₂ O)	Oligo (n)
25	2.5 μL	10 μM	20.0 μL	200 μM	2.0 μL	0.5 μL	200 pmol
25	2.5 μL	50 μM	20.0 μL	1 mM	2.0 μL	0.5 μL	1.0 nmol
25	2.5 μL	100 μM	20.0 μL	2 mM	2.0 μL	0.5 μL	2.0 nmol
100	10.0 μL	10 μM	80.0 μL	1 mM	1.6 μL	8.4 μL	0.8 nmol
100	10.0 μL	50 μM	80.0 μL	2 mM	4.0 μL	6.0 μL	4.0 nmol
100	10.0 μL	100 μM	80.0 μL	10 mM	1.6 μL	8.4 μL	8.0 nmol